

**CPAP Nasal Prongs:**

- ❖ Available in seven different sizes based on the patient's nares. Select the appropriate prong size by considering the baby's weight and nares size.
- ❖ To determine the correct prong size, position the dots on the sizing guide over the infant's nostrils.

<b>CPAP Nasal Prong size</b>	<b>CPAP Nasal Prong Number</b>	<b>Baby Weight (g)</b>	<b>Component Weight (g)</b>
Micro	98-00-116	Up to 500	1.1
Small	98-00-115	500-1000	1.1
Medium	98-00-110	1000-1500	1.1
Medium Wide	98-00-111	1000-1500	1.2
Large	98-00-112	1500-2200	1.2
Large Wide	98-00-113	1500-2200	1.4
X-Large	98-00-114	2200-3000	1.1