

## Eye Pad

- ❖ Select the eye pad size based on the baby's head circumference and eye coverage. The eye pad must fully cover the eyes without exerting pressure.
- ❖ Adjust the elastic/hookable strap around the head

### a) Elastic

Size	Part Number	Baby Weight
Small	76-00-030	1000 – 1500 grams
Medium	76-00-031	1500 – 2000 grams
Large	76-00-032	> 2000 grams

### b) Hookable

Size	Part Number	Baby Weight
Small	76-00-059	1000 – 1500 grams
Medium	76-00-060	1500 – 2000 grams
Large	76-00-061	> 2000 grams