

Bubble CPAP Head bonnet & Measurement Tape

- ❖ Eight types of head bonnets are available. Measure the baby's head circumference using the Head Measurement Tape (Part No. 73-00-037) to select the appropriate head bonnet.
- ❖ Use the hookable lace for looping and attaching the patient interface with head bonnet.

Head Bonnet size	Length	Part number
XX Small	17 to 19 cm	76-00-003
X Small	19 to 21 cm	76-00-006
Small	21 to 23 cm	76-00-007
Medium	23 to 25.5 cm	76-00-008
Large	25.5 to 28 cm	76-00-009
X Large	28 to 30 cm	76-00-010
XX Large	30 to 33 cm	76-00-011
XXX Large	33 to 36 cm	76-00-012