

MAXFLOW KIDS NASAL CANNULA

Part Number	50-05-261 XS (Blue)	50-05-113 S (Red)	50-05-114 M (Yellow)	50-05-115 L (Purple)	50-05-116 XL (Green)
Patient type	Premature	Premature	Neonate	Infant	Paediatric
Flow Range	4 to 8 LPM	4 to 9 LPM	4 to 10 LPM	4 to 15 LPM	4 to 15 LPM
Wangpad Sizes	76-00-057		76-00-058		
Expected to fit patient Weight in kg)	1-2	1-3.5	2-8	4-15	7-25
May fit patient Weight in kg)	0.5-1 & 2-2.5	0.8-1 & 3.5-4	1-2 & 8-10	3-4 & 15-20	5-7 & 25-30

HIFLEX KIDS NASAL CANNULA

Part Number	50-05-268 XS	50-05-123 S	50-05-124 M	50-05-125 L	50-05-126 XL
Patient type	Premature	Premature	Neonate	Infant	Paediatric
Flow Range	4 to 8 LPM	4 to 9 LPM	4 to 10 LPM	4 to 15 LPM	4 to 15 LPM
Wangpad Sizes	76-00-057		76-00-058		
Expected to fit patient (Weight in kg)	1-2	1-3.5	2-8	4-15	7-25
May fit patient (Weight in kg)	0.5-1 & 2-2.5	0.8-1 & 3.5-4	1-2 & 8-10	3-4 & 15-20	5-7 & 25-30

- The Maxflow Kids and Hiflex Kids cannula are available in five sizes based on patient type and weight; select the appropriate size accordingly.

MAXFLOW ADULT NASAL CANNULA

The Maxflow Adult nasal cannula is available in three sizes based on the flow range; select the appropriate size accordingly.

Performance Specifications	50-05-195	50-05-196	50-05-197
Size	Small	Medium	Large
Flow range	10 to 40 LPM	10 to 60 LPM	10 to 60 LPM
Insp. Resistance at 60 LPM	> 13 cmH ₂ O		
Leakage at 60 LPM	> 9 LPM		