

T-bar Prongs:

- ❖ Available in Five different sizes based on the patient's nares and baby weight Select the appropriate T-bar prong size.

CPAP Nasal T bar Prong size	CPAP Nasal T bar Prong Number	Baby Weight	Component Weight
Size 0	98-00-129	<700 grams	2.3 grams
Size 1	98-00-130	<700 grams	2.3 grams
Size 2	98-00-131	700-1250 grams	2.3 grams
Size 3	98-00-132	1250-2000 grams	2.7 grams
Size 4	98-00-143	2000-3000 grams	2.8 grams